

Model: ISO-Poultry1 Isolator

Heavy duty flexible film Isolator to house chickens available to operate both in positive and negative pressure

Dimensions

Overall (height adjustable): 110-210cm(L) x 85-110cm(W) x 120-160cm(H)

Containment Unit: 100-200cm(L) x 95cm(W) x 745cm(H)

Port options:

- Standard Quadro-Lock port 42cm with two easy open doors which provides extra security
- Standard Port 42cm with flexible film covers
- XL Quadro-Lock port 46cm with two easy open doors which provides extra security
- XL Port 46cm with flexible film covers

Canopy:

Heavy duty PVC for fitted with soft detachable sleeves for user comfort. To avoid access of animals into the sleeves, small lids are fitted. Position of sleeve will be manufactured according to request.

Additional options - like hooks and pockets - inside and out are optional.

Air flow:

Fan unit with tubing and valves to allow controlled airflow independently for each Isolator to work as a self-sufficient unit.

Filters:

Air –inlet and - outlet with fitted HEPA Filters for highest security even if airflow would be interrupted.

External Frame:

Frames are available in stainless steel or enamelled aluminium frame. Available as fixed height or adjustable frames to suit your needs

Accessories:

NKP is happy to provide accessories like food hopper, drinking container and enrichment as you require.

Delivery: Since each unit is custom manufactured to your specification, please allow 9-11 weeks for manufacture - delivery time according to location

Assembly: We offer delivery and assembly of the units in our service package. Whilst there are detailed instructions are available we do recommend that assembly is performed by trained personnel.

Isolator Specification Sheet

Dimensions

Overall (height adjustable): 110-210cm(L) x 85-110cm(W) x 120-160cm(H)

Containment Unit: 100-200cm(L) x 95cm(W) x 745cm(H)



Special feature:

sleeve lids to restrict access of the animals from the inside.

Optional hooks and pocket can be added to your canopy.



Large Port option

this picture shows the half/half XL-Quadro-Lock port option with a solid door on the inside and flexible cap on the outside

