

## Model: Hydraulic double tier

This hydraulic Isolator combines mobility with ergonomic working heights at a very compact and space saving design.

Different Port and Canopy options are available. This Isolator is ideal for experimental set-up or small colony breeding.

The overall dimensions can be changed to fit your needs.

### Specifications:

- Overall height 165cm up to 230cm
- Hydraulic lift up to 70cm
- Canopy length between 1m and 1.5m
- Canopy height 70cm - 85cm (overall height requirements need to be taken into account)
- Frame depth 80cm - 100cm
- Frame on castors for mobility
- Hepa filtration for Inlet and exhaust
- Operation can be switched between positive and negative pressure
- Castor base with legroom to allow working with a chair

